
移动高压GMP清洗机 MOBILE HIGH PRESSURE GMP WASHER

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上海弋凌流体机械设备有限公司 Shanghai Yiling fluid machinery equipment Co., Ltd

INNOVATIVE SOLUTIONS

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Cleaning unit pharmaceutical grade, suitable for washing bins, drumsand bio reactors (with or without agitators). It provides reliable, repeatable and verifiable results to meet the requirements of the pharmaceutical industry.

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The hydrokinetic nozzle provides a 360° movement over a defined time period. It is automatic with a water consumption of < 40L/min. The cleaning time varies between 5 to 10 minutes for full coverage. The water and time reduction for cleaning a bioreactor will be reduced by >90% compared with manual cleaning

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LAYOUT OF A SIMPLE CLEANING STATION

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
TECHNICAL SPECIFICATIONS

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Structure 结构

Frame and housing made of stainless steel AISI 304, with noise and heat reduction, equipped with 4 wheels.

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High pressure pump designed for the direct connection to the facility utilities, pump head made of stainless steel AISI 316.

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Washing system (nozzles) 洗涤系统(喷嘴)

2 connections to the hydrokinetic nozzles and/or lances, specifically designed for the application.

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Compressed air connection 压缩空气连接

Filtered process air for blowing out the piping system 过滤工艺空气,用于吹出管道系统

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iemens SPS S7 with integrated HMI (Touch Screen). Graphical visualization of the de process and the process parameters.

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Water connection (customer) 水需求 (用户)

Pressure 1.5–3 bar Flow max. 40 l/min Temperature max. 90°C

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
移动高压GMP清洗机 MOBILE HIGH PRESSURE GMP WASHER


CLEANING PATTERN

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COMMISSIONING, IQ/OQ AND FAT/SAT 调试、IQ/OQ和FAT/SAT

The cleaning machines are completed in the factory, installed on the in-house test bench, started up and calibrated by Eling specialists. As part of the internal factory test (pre-commissioning), the conformity of the machine is demonstrated in documented way. During a factory acceptance test (FAT), the cleaning machine can be checked and accepted. The commissioning on site as well as the training of the operating and maintenance personnel is carried out by qualified Eling specialists. Installation and Operational Qualification (IQ / OQ) and also Performance Qualification (PQ) support can be provided based on defined test specifications.

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The machine is controlled via a modern Siemens Control-system which stands for reliability, flexibility and worldwide support. The visualization (HMI) runs on a Siemens TP Comfort Panel. The cleaning cycles are individually configured to be adapted to the soiling and the specific requirements of the wash-load.


The integrated user management on the software allows easy configuration of the access-levels of each user-group. Specific user-group have different access-rights; therefore, configuration or changes on cleaningprocesses can only be made by authorized personal.

Process data monitoring and trending of relevant parameters is available in the control concept and can be set individually, depending on the washing recipe.

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INTEGRATION IN CUSTOMER NETWORKS

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Integration into customer networks can be executed in different ways. Possible interfaces are Siemens Profinet, Ethernet and USB. MES System-connection (manufacturing execution system), via acustomer-side DCS System (distributed control system), eg for sending or receiving the cleaning recipes. MES System-connection via SCADA-functionality (Supervisory control and data acquisition) by Siemens software. Document-network for data-storage on a fileserver with active directory for machine-generated pdf. PAN-connection (Personal area network), for time-synchronisation, user-log-in and machine generated protocols; OPC (OLE for process control).

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Profinet-Communication

让清洗符合规范要求! Make the cleaning meet the requirements of the specifications!

上海弋凌流体机械设备有限公司 Shanghai Yiling fluid machinery equipment Co., Ltd.

Tel:+86-021-57480503 Web:www.yl-fluid.cn 上海市奉贤区航南公路 5088 弄9-1501 9-1501, Lane 5088, Hangnan Road, Fengxian District, Shanghai